

Fig 2 Illustration of Spark-DAT principles




Fig 3 Comparison of clean and dirty steel samples




Fig 4 Example of inclusion counting with the algorithms Peaks and Composition

concentration of 0.075% is observed. This is illustrated in Figure 2 for Al₂O₃ inclusions in an Fe matrix.

The intensity of the baseline signal is proportional to the concentration of Al atoms dissolved in the matrix and the intensity of a peak depends on the amount of Al atoms contained in the inclusion(s) ablated by a single spark. As a consequence, the number of peaks is related to the number of such inclusions and their intensity to factors such as the inclusion size and the concentration of Al in the inclusion.

PRACTICAL DETAILS AND BENEFITS

Inclusion analysis can be performed stand-alone in about 7secs for a single measurement. However, combined with elemental concentration analysis, it offers more possibilities, while increasing the time necessary for a single measurement to about 22secs only (against 19secs for the standard OES analysis).

The main benefits of such a short combined analysis are:

- Inclusion analysis during steel elaboration on most of the samples analysed normally by OES
- Quantitative inclusion size determinations
- Quantitative determination of the total oxygen content in killed steels

The standard, quick OES sample preparation can be used for Spark-DAT inclusion analysis. However, with paper grinding, the paper should be chosen in order to avoid any influence on the analysis of the inclusions of interest (eg, using SiC paper when Al₂O₃ inclusions have to be analysed). For advanced Spark-DAT, milling is necessary to guarantee quantitative results.

Other benefits of using Spark-DAT inclusion analysis are:

- Low investment costs, the spectrometer being sufficient for performing inclusion analysis in addition to analysis of elemental concentrations
- No additional cost and time for operations compared to the standard OE spectrometer. Maintenance, service and sample preparation for inclusion analysis are equivalent

EVALUATION OF NUMBER AND TYPE OF INCLUSIONS

The algorithm Peaks allows evaluating the number of inclusions counting intensity peaks on the channels of element present in inclusions. A peak is defined as an intensity signal I_{peak} higher than a threshold situated at the mean intensity m of the element dissolved in the matrix plus three times its standard deviation SD :

$$I_{\text{peak}} > m + 3 \cdot SD_{\text{matrix}}$$

As shown in *Figure 3*, steel samples can easily be classified as clean or dirty according to the number of peaks counted on the channels of the inclusion elements.

The algorithm *Composition* allows counting coincidences of peaks, ie, peaks appearing on the channels of several elements consecutively to the same single spark, as illustrated in *Figure 4*. The coincidence of a peak on Ca and S channels signifies that these two elements are part of the same inclusion, very often a calcium sulfide (CaS) inclusion. Coincidences of up to four channels can be counted with this algorithm. This enables the chemical formulation of complex inclusions or inclusion clusters. Furthermore, the possibility to check for non-coincidences on top of coincidences helps removing ambiguities between inclusions that share elements, like Al_2O_3 and $\text{Al}_2\text{O}_3\text{-CaO}$ or CaO and $\text{Al}_2\text{O}_3\text{-CaO}$ as shown in *Figure 4*.

QUALITATIVE INCLUSION SIZE ANALYSIS

Knowing the size of the inclusions and in particular their size distribution is important because large inclusions are normally the most detrimental to metal quality. The algorithms *Peaks* and *Composition* can be used to count inclusions (signals) belonging to different size (intensity) classes. Setting the threshold 3-SD above the intensity of the element in the matrix allows us to count all the inclusions that are large enough to be detected. Setting it higher, for example at 9 or 15-SD, as in *Figure 5*, allows us to count inclusions with a larger size.

Calculating the inclusions between consecutive threshold values gives the number of inclusions in their specific size class. In the example, the number of peaks and coincidences between 3 and 9-SD correspond to small size inclusions, between 9 and 15-SD to medium size inclusions and higher than 15 to large size inclusions. Such calculations are possible with more than three size classes. This allows generating inclusion size distribution diagrams that are qualitative, but can be extremely useful for example in order to compare the inclusion distributions in samples of different heats.

QUANTITATIVE INCLUSION SIZE ANALYSIS

Nowadays, measurements of inclusion size and size distribution are normally performed by techniques such as optical microscopy and scanning electron microscopy (SEM/EDX). However, the methods are time-consuming, taking typically 2-4hrs per analysis in a very competent laboratory.

The algorithm QuIC (Quantification of Inclusion Content) enables quantitative analysis of inclusions in term of size and size distribution. With it, average ESD of inclusions belonging to several size classes can be calculated.

Figure 6 presents an example of size distribution




Fig 5 Example of inclusion counting using different threshold values

diagram, based on the application of the QuIC algorithm and size calculation method for several production samples.

Other parameters, such as the surface fraction S_f (see *Figure 7*), of a given inclusion type and the total oxygen content, can also be calculated from the results of the algorithm QuIC.

QUANTITATIVE ANALYSIS OF OXYGEN

Due to the relatively low sensitivity of the analytical line, the standard routine quantitative analysis of total oxygen by OES is currently semi-quantitative below 60ppm and very limited below 30ppm. In the steel industry, oxygen at very low concentration is therefore normally analysed by dedicated combustion analysers.

For killed steels treated with the typical deoxidisers (Al, FeSi, FeMn and Ca, etc.) this type of size determination is quantitative, because most of the oxygen is insoluble, more precisely, in the form of oxide inclusions. With Spark-DAT the total oxygen concentration can therefore be calculated directly from the amount of oxygen contained in the oxide inclusions obtained with QuIC.

The method is quantitative and determines the oxygen concentration by OES well below 30ppm, as demonstrated in *Figure 8*, which compares results obtained with those of a combustion analyser for low alloy steel samples taken from the casting mold (*samples and combustion results with permission of R Dumarey and F Medina, ArcelorMittal, Gent, Belgium*).




Fig 6 Example of size distribution of Al_2O_3 for three alloy steel samples




Fig 7 Example screen display




Fig 8 Comparison of Spark-DAT and combustion analysis results for samples coming from the casting mold

INCLUSIONS QUANTIFIABLE WITH SPARK-DAT

Various endogenous and exogenous inclusions can be observed directly or indirectly in steel using the ARL 4460 with Spark-DAT, eg, oxides (Al_2O_3 , MgO , CaO , MnO , TiO_2 , SiO_2), spinels (Al_2O_3 - CaO , Al_2O_3 - MgO) and sulfides (CaS , MnS , AlS). Detection of an inclusion type is limited mainly by the sensitivity of the analytical lines used, the size of the inclusion and the concentration of the inclusion elements as soluble elements in the matrix: higher line sensitivity and lower soluble content determine smaller inclusion ESDs. For example, for a steel containing 50ppm of Al the minimum detectable ESD is about $1\mu m$, while for 0.2% of Al the minimum is $4.5\mu m$.

CONCLUSIONS

The Spark-DAT option of the ARL 4460 metals analyzer increases the versatility of the OE spectrometer. From routine use to research, Spark-DAT provides quick, simple and cost-effective solutions for inclusion analysis in the steel industry. Among all the inclusion analysis techniques available today, the Spark-DAT based methods are the fastest. It allows ultra-fast online counting and identification of inclusions, determining their size and total oxygen content, ranging from several seconds to a few minutes, making it highly effective for controlling inclusions and steel cleanliness during production. The sample and its surface preparation, as well as instrument maintenance and consumables, are equivalent to a standard OES instrument, milling being required for quantitative analysis. This signifies extremely low operating costs compared to other techniques that need dedicated instruments. Also, the ability to obtain elemental analysis information and inclusion contents with a single OES drastically reduces investment and operating costs.

ACKNOWLEDGEMENTS

We are indebted to Rolf Didriksson from SSAB (Oxelösund) for enlightening (sparking) discussions. We acknowledge Ronny Dumarey and Frank Medina from ArcelorMittal (Gent) for extremely valuable help and discussions, and for giving permission to publish results obtained on their samples. **MS**

REFERENCE

[1] J-M Böhlen & R Yellepeddi, "Combined quantitative analysis and ultra-fast analysis of non-metallic inclusions by optical emission spectrometry", *Millennium Steel* 2009, pp167-171

Kaizhen Li is Application Specialist, Jean-Marc Böhlen is Product Manager and Edmund Halász is Manager R&D, all for OES products at Thermo Fisher Scientific, Ecublens, Switzerland.

CONTACT: analyze@thermo.com